


THE SUMMIT WS-19 Massage Chair Owner's Manual



www.wellnessshop.ca/massage-chairs/

Please read SAFETY INSTRUCTIONS carefully before using the chair.

**The Summit
WS-19
Massage Chair**

Introduction

1. The Summit's intelligent 3D massage helps you enjoy the feeling of having a professional masseuse right at home.
2. Using its automatic ShapeDetect system at the beginning of each session, the Summit delivers an optimized massage for people of all sizes.
3. With its heating massage capabilities, the Summit can relieve pain and improve blood circulation.
4. The Summit's LongRail technology provides continuous therapy with its massaging rollers starting from the head, down the back, around the glutes, and close to the back of the knees.
5. Easily place the Summit in ZeroG mode to evenly disperse your weight throughout your spine and alleviate back pain.
6. Extend the legrest and adjust the footrest to any angle you need.
7. The backrest and seat frame move together seamlessly, allowing you to take full advantage of surrounding space.
8. Plug in your mobile phone via the chair's USB port to listen to music while you enjoy the massage.
9. The Summit's special air pressure design will promote blood circulation and restore vitality and energy.
10. The remote control will display various massage parts, methods, and status updates to help you optimize your time in the chair.

- ✦ 1. Thank you for choosing the Summit WS-19 Massage Chair.
 - ✦ 2. Please read this manual carefully before using the chair.
 - ✦ 3. Please keep this manual in a safe place for future reference.
 - ✦ 4. Please pay attention to the SAFETY INSTRUCTIONS.
- ✦ Note: We reserve the right to design modification, which is subject to change without notice. The pictures in this manual are for reference only. Please refer to the actual product.

Preparation	Safety instructions	1
	External structure parts	6
	Internal structure parts	6
	Mechanical components diagram	7
	Mechanical components list	6
Installation	Installation	9
Massage Operation	Location placement	17
	Remote control	19
	Operation instructions	23
	After massage	25
Maintenance	Cleaning	26
Others	Troubleshooting	27
	Technical data	28

Please keep this manual in a safe place for future reference.



Safety instructions

- Please read this manual carefully to make sure the installation and operation of the massage chair is correct.
- The "Safety Instructions" section contains both warning and cautions.

Please pay careful attention to the following symbols:

	Warning	Possible death or injury could occur
	Caution	Possible injury or object damage could occur

Symbol meaning

	Prohibited behavior		Prohibited disassembly
	Follow these instructions to operate		Correct operation of power cord

Keep these instructions for future reference.

Preparation

Safety instructions



Warning


	<ul style="list-style-type: none"> • People with osteoporosis and/or heart disease, including pacemakers or other electronic medical devices. • People with high fever or who feel unwell. • People with certain perceptual disabilities. • People over 80 years old, and children. • People who are sensitive or allergic to heat. • People who are intoxicated. • Pregnant women should consult their doctor before use. • Do not sleep in the chair. 	People who must not use the chair
	<ul style="list-style-type: none"> • Make sure no children, pets, or objects are under, behind, or in front of the chair when reclining or retracting/extending the legrests. • Do not allow children or pets to play on the chair, especially during operation. This includes standing or sitting on the back or armrest. • If your body feels abnormal during the massage, stop immediately and consult your physician. • Please do not use a strong massage setting for your head. • Do not massage your belly or knees, and be careful when massaging your neck. • Do not use when your body is wet. 	Accident/injury warning
	<ul style="list-style-type: none"> • Be sure to insert the power cord completely into the correct power socket to reduce the chance of short circuiting and starting a fire. • Unplug the chair when not in use. Stop all massage functions and press the OFF key before unplugging. 	Fire/electrical shock warning
	<ul style="list-style-type: none"> • Do not use the chair when the ventilation is blocked. • Do not fold the chair's cushions for storage. • Do not operate the chair when no one is sitting in it. 	Accident warning
	<ul style="list-style-type: none"> • Do not use if the power socket is broken or loose. • Do not use with any other voltage except the rated voltage. • Do not use if the power cord has been damaged, twisted, or knotted. 	Fire/electrical shock warning
	<ul style="list-style-type: none"> • If the chair is broken, dysfunctional, or has been exposed to water, call a professional repairman or send it to an assigned service agent. Do not modify, disassemble, or do your own maintenance. • If the power cord is damaged, do not replace it yourself. Contact your service agent for a replacement. • To avoid electric shock or injury, do not open or disassemble any part of the chair. • Make sure no parts of the leather or cloth material are damaged. If you come across even a small tear, stop using the chair, unplug it, and have it repaired. 	Disassembly and maintenance

Preparation



Preparation

 Caution	
<ul style="list-style-type: none"> Be sure to set the chair on level ground. Position the back of the chair 5' from the wall to allow for sufficient space to recline up and down. 	To ensure a full recline and avoid falling suddenly
<ul style="list-style-type: none"> The suggested operating time is 20 minutes total. Spend 5 minutes on each position, beginning with slight-step massage. 	To avoid excessive massage, which may cause soreness
<ul style="list-style-type: none"> Before sitting on the chair, make sure the massage heads are in the correct position. Make sure to place the remote control in the pocket or on the bracket. 	To avoid accidents or injury
<ul style="list-style-type: none"> Check that the power switch is on the OFF position before plugging in the power cord. 	To avoid accidents or fire hazard
<ul style="list-style-type: none"> Turn off all controls before removing plug from outlet. Make sure plug is completely inserted into the outlet. 	To avoid accidents or fire hazard
<ul style="list-style-type: none"> When checking for wire damage, beware of water or debris on or near the wire. Make sure the plug is inserted into a stable socket that is not bent or cut off.  	To avoid electric shock
<ul style="list-style-type: none"> Don't press hard against the massage heads, or put your hands and feet between them, during operation. Don't put hands or feet in the gap of the chair's mechanical parts. Don't drop anything into the chair slot. Don't stand on, sit on, or set anything heavy on the armrest, legrest, and backrest, especially when the chair is reclined. Don't move the chair while in operation. Don't put your hands or head underneath the legrest. Don't wear hair accessories while using the chair. Don't expose bare skin, or wear thin clothes, as this may cause irritation or cryogenic burns. 	To avoid accidents or injury
<ul style="list-style-type: none"> Keep the chair away from sharp objects, as well as wet, dusty, corrosive, inflammable, or explosive environments, including stoves and other heating appliances. Never use an electric blanket or other heating/ electrical product while using the chair. Never place the chair in a small space with bad air circulation. Never place the chair outdoors or other high-temperature environments. 	To avoid fading or discoloration to upholstery or hardened deterioration on PU

 Caution	
<ul style="list-style-type: none"> Don't use or store the chair in a wet environment, such as near a swimming pool or bathroom. Keep the chair and the remote control dry at all times. 	To avoid electric shock
<ul style="list-style-type: none"> Make sure the ground wire is installed correctly. If you don't know how to install it, hire a professional electrician. Do not use a plug adapter. If the plug does not fit into the socket properly, do not change the plug by yourself. Hire a professional electrician. 	To avoid electric shock
<ul style="list-style-type: none"> Unplug the power cord before cleaning the chair. Unplug the power cord after turning off the chair. Unplug the power cord when there's any damage. 	To avoid electric shock or fire hazard
<ul style="list-style-type: none"> Try not to pull the wire when unplugging the power cord. 	To avoid accidents
<ul style="list-style-type: none"> Turn off the power switch and unplug the power after using the chair. 	To avoid damage or accidents
<ul style="list-style-type: none"> Unplug the chair immediately after losing electricity. 	To avoid damage or accidents

Preparation


Grounding Information

✦ The chair is equipped with a power cord that has a grounding conductor. It must be plugged into a socket that fits appropriately and is properly installed.

Fuse Replacement

✦ When the fuse burns out, turn off the switch and unplug the power. Open the lid covering the fuse, located in the switch box. Replace with the same size fuse, then cover the lid. Or, hire a professional electrician.



 Ground wire is prohibited in these places	<p>Gas pipe: Can cause explosion or fire</p> <p>Phone wire or lightning rod: Can cause electric shock or fire during thunderstorms</p> <p>Water pipe: Ground wire is useless with plastics</p>
---	---




5 External structure parts


Internal structure parts

6

Preparation




Preparation



7

Mechanical components diagram



Preparation

Mechanical components list

8

No.	Name	Qty	No.	Name	Qty
1	Remote control	1	22	Pillow	1
2	Right side fender	1	23	Right shoulder massage device	1
3	Air pump	2	24	Right armrest	1
4	Connecting rod	1	25	Seat frame surface	1
5	Legrest & footrest	1	26	Right side of seat airbag	1
6	Front fender	1	27	Seat support bar	2
7	Left side fender	1	28	Castor	2
8	Left side seat airbag	1	29	Power case	1
9	Backrest electric actuator	1	30	Right armrest magnetic valve	1
10	Bottom seat frame	1			
11	Footrest electric actuator	1			
12	Sliding rail	2			
13	Top seat frame	1			
14	Left armrest	1			
15	Electric magnetic valve	2			
16	Left shoulder massage device	1			
17	Massage machine	1			
18	Back cover	1			
19	Back frame	1			
20	Backrest front cover and seat plate	1			
21	Back and seat cushion	1			


Preparation



Step 1. Remove chair and accessories from the box


- ① Make sure that two people are removing the accessories and chair.
- ② Confirm all accessories are there.

Hold the metal seat frame and take the chair out of the box.



Message chair main body

The accessories



Head pad

Back seat cushion

Remote control bag



Right side airbags of the seat (zipper side upwards)

Left side airbags of the seat (zipper side upwards)



Remote control



Power cord



Thigh shiatsu panel



Left shoulder massage device



Right shoulder massage device

- ③ Take out both armrests and place them gently on a flat surface.



Right armrest

Left armrest

- ④ When taking out the legrest device, hold both the front and back two sides.




Legrest device

Installation

● Check the surrounding environment for obstructions

At least 12 inches from the wall At least 4 inches from the wall



Be sure to place the chair on a blanket or other soft surface area to avoid ruining your floor due to its heavy weight.


Caution:

- Don't place the chair near wet areas, such as a bathroom.
- Place the chair on a flat surface.
- Don't place the chair under direct sunlight or near any heating source, such as a stove.

Installation

Step 2. Return massage chair to default position


- ① Insert the power cord into the chair and then the wall outlet.
- ② Set power switch on the chair to ON.
- ③ After the chair returns to its default position, set the power switch to OFF and unplug the chair from the wall outlet.



Step 3. Assemble the legrest

- ① Unscrew the bolt from the hanger of the legrest. ② Connect the three air pipes on the legrest with the air tap on the main body according to the colors. ③ Connect the power cord to the power socket.
- ④ Raise the rotation shaft board on hangers, raise the legrest, then put the rotation shaft into the hangers of the main body, and place the shaft board back down. ⑤ Screw the bolt back in tightly.

⚠ Caution: Make sure the air tubes and wires do not get twisted.




Installation

Step 4. Assemble the left and right armrest

- ① Unscrew the M6 bolts on the seat frame. ② Press the button (see below) inside the back of the armrest, then insert the bolt of the side backrest. ③ Insert the square tube into the side of the seat frame. ④ Match up the front and back panel bolt holes on the metal seat frame, and screw the m6 in. ⑤ Connect the air pipe of the same color onto the armrest.

⚠ Avoid squeezing the air pipe during assembly

✦ To disassemble the armrest: Disconnect the air pipes between the seat and armrest, undo the four screws on the seat, lift up the front of the armrest, then let the armrest square tube move out from the seat frame plate. Press the button (see below) to pull the end of the armrest and at the same time, move the bolt hole away from the bolt on the side of the backrest.




Installation



13 Installation

Step 5. Install the left and right seat airbags and the remote control

- ① Align the two cylindrical pins to the seat frame holes, then press down and make sure the connection is tight.
- ② Connect the air pipe on the airbag with the air tap behind the seat.
- ③ Attach the remote control bag to the chair, between the right armrest and airbag.
- ④ Insert the remote control into the right armrest pocket, and plug in the cord.




Installation

Installation 14

Step 6. Connect the cushions

Connect the back seat cushion and head pad with the zipper and Velcro.
Make sure everything fits correctly in order for the body tests to be accurate.




Installation

Step 7. Install left and right shoulder massage device

▲ Be sure to assemble the left shoulder massage device with the left side of the backrest, and right shoulder massage device with the right side of the backrest.

- ① Unscrew the bolt affixed to the inside of the back armrest.
- ② Connect the speaker wire and shoulder massage device's air pipe with the side of the backrest.
- ③ Align the two axle holes of the shoulder massage device with the two axles on the backrest, then connect them.


▲ During installation, do not squeeze the air pipe and speaker wire to avoid damage.



15 Installation

Step 7. Install left and right shoulder massage device

- ④ Turn the shoulder massage device downward, until it no longer moves.
- ⑤ Unzip the airbag of the shoulder massage device, then lift up the airbag and screw the device into the armrest.
- ⑥ Repeat these steps for the other shoulder massage device.



Installation


Installation 16

Step 8. After assembly

- ① Plug in the power cord and remote control cord.
- ② Turn on the power box with the O/I button.
- ③ Use the USB connection to play music, if desired.
- ④ Use the remote control to start the massage.



Make sure the power switch is off before you plug the power cord into the wall socket.



Installation



Caution


Don't place the chair near wet places, such as bathrooms and swimming pools.

Don't place the chair under direct sunlight, or near other heat sources such as stoves.

Place the chair on an even surface.

Check area around the chair

Ensure there are no people, pets, or objects around the chair before turning it on.



Check power cord

Be sure the power cord isn't dusty, bent, or damaged.



Measure the surrounding area


Make sure the chair is at least 4 inches from the wall.



Make sure the front of the chair has at least one foot of clearance.

Check the surrounding area

When extending the backrest or legrest, be sure to check the area is clear.



Massage operation

Caution

Don't move chair when it's working to avoid accidents or damage.

How to move

- Bind the cords and place them in the seat to avoid damaging them when moving the chair.
- If you plump the chair, it might damage the inner parts of the chair.
- Don't use the castor if the floor is wooden or delicate in any way.
- Make sure the power is off and the power cord doesn't touch the floor.



Lift up the legrest and push the chair to move it to its destination.



Move the chair as shown in the picture to avoid any damage.

Before adjustment

Make sure the rollers are under the seat before sitting down.



If the acupuncture detection test does not detect your head, the chair will default to the middle position.

Correct position: Your back should be snug to the chair back, with your head snug to the pillow.


The remote control will beep when the acupuncture detection test is completed. Press UP/DOWN when needed.




Suggestions

If the massage is not intense enough, take off the pillow or cushion.

If the leg massage is not strong enough, place the Shiatsu plates in the leg area as shown.




Because the chair is heavy, it may ruin your floor if it stays put for a long time. Place a blanket or other soft items under the chair to protect your floor.



















Massage operation



Massage operation



-  When the chair is on standby, press ON/OFF to turn it on. Then choose auto or manual massage. The chair will run a shoulder height detection test and then begin the massage. Press ON/OFF again to turn the chair off, and the chair will return to its default storage position.
-  Press this button to pause the current action, and press it again to continue.
-  Press this button to switch to auto massage mode. Choose from Recovery, Pain Relief, Stretch, Music Sync, and Experience.
-  Press this button in the main interface to switch to the menu, where you can choose different massage options.
-  Press this button to turn your music on or off or take a phone call.
-  Use this to navigate the menus or adjust the volume of your music.
-  Press this button to confirm menu choice and enter the next menu.
-  Press this button to adjust width (wide, medium, and narrow.) This can only be adjusted in manual mode, except when using the kneading or kneading & knocking technique.
-  Press this button to adjust the speed between level one, two, and three. This can only be adjusted if in manual mode, except when using the finger press technique.
-  Press this button to adjust the intensity of the air massage between level one, two, and three. This can be adjusted in auto or air massage mode.
-  Keep pressing this button to make the legrest go up. Release it to lock it in position.
-  Press this button to adjust the speed of the foot rollers between fast, slow, or off. You can also turn on the calf massage or turn it off.
-  Keep pressing this button to make the backrest go up and the legrest go down. Release it to lock it in position.
-  Keep pressing this button to make the legrest go down. Release it to lock it in position.
-  Press this button to put the chair in ZeroG position. Pressing it again will restore the chair to its default position.
-  Keep pressing this button to make the backrest go down and the legrest go up. Release it to lock it in position.

Massage operation



Main interface

After turning on the remote control, the LCD display will appear after a few seconds, showing the massage position, methods, duration, etc. This information will be updated in real time.

Massage function

- Kneading
- Tapping
- Knocking
- Pressing
- Zero gravity
- Massage time left

Heating

Music

Message position

- Message speed (slow, mid, fast)
- Message width (narrow, mid, wide)
- Air intensity (weak, mid, strong)
- 3D strength (4 levels)

Menu interface

Auto massage: Press the Menu button, then press OK to enter automatic mode: Recovery, Pain Relief, Stretch, Music Sync, Experience, and Special. The Special mode includes five additional automatic massage programs.

- White-collar
- Shopaholic
- Back Yoga
- Athlete mode
- Care of aged

Manual message: Press the Menu button, select manual message and press OK. Choose among Mode, Position, Width Adjustment, Speed Adjustment, and 3D, and press OK to select the desired functions in each.

When the chair is in manual fixed point mode, use the UP/DOWN button to adjust the position.

Note:

- When selecting a massage method, the default position is the whole body.
- The part and fixed point massage do not need shoulder height detection. When you switch to a whole body massage or similar, the shoulder height detection test will begin.

Manual message menu:

- Auto
- Manual
- Air
- Roller
- Heating
- Settings

Manual message sub-menu:

- Mode
- Position
- Width
- Speed
- 3D

Manual message sub-sub-menu:

- Overall
- Partial
- Foot part
- Knock
- Knock&Tap
- Shiatsu

Air massage: Press the Menu button, select air pressure, then press OK. Select Position or Intensity to adjust the function accordingly and then press OK.

Roller massage: Press the Menu button, select roller, and choose from Quick, Soft, and Stop.

Note:

- For most types of massage, the foot roller will be on by default. If you don't want it, press the Stop button in the foot roller menu.
- For the hand & shoulders and waist & seat air pressure massages, you will need to go into the foot roller menu to turn it on.

Heating: Press the Menu button, select heating, and choose from various heating functions. The heating signal will flash on the remote control.

Settings: Press the Menu button, select settings, and choose Timer, Language, or Stop.

Timer: Choose 10, 20, or 30 minutes. When the time is up, the chair will power down automatically.

Language: Choose from English, Chinese, or Korean.

Stop: Press this button to stop all massage function.

The first time using an auto program, the chair will start acupuncture position detection, which will be indicated on the remote control screen. After it is finished, the remote control will buzz five times. Use the UP/DOWN button to readjust the chair position if needed.

Acupuncture point tracking

Massage operation


Massage operation



23 Operation instructions

Operation instructions 24

Massage operation




Step 1. Turn on power

- ① Plug end of power cord into A.
- ② Plug remote control cord into B.
- ③ Plug power cord into wall socket.
- ④ Use the power switch to turn on the chair.
- ⑤ Press ON/OFF on the remote control to turn it on.

Step 2. Massage

- ① Follow the instructions on the remote control.
- ② When the massage ends, "shutting down" will be displayed on the remote as the chair returns to its default position.



Step 3. Turn off the power

- ① Press ON/OFF to turn off the chair, and it will restore automatically.
- ② After the rollers return to storage position, turn off the power switch on the bottom of the back of the chair.
- ③ Unplug the power cord from the wall socket.

Warning
Before using the chair, be sure to check for any damage or deterioration. Only use a proper power outlet.

Attention
1. Before sitting, make sure nothing is pinched between the armrest, legrest, backrest, and chair seat.
2. Make sure the chair is restored to its default position before sitting.
3. Do not stand on the chair.

Massage operation

● Legrest up

Press this button on the remote control to make the legrest go up. Release the button to lock in position.



● Legrest extend and retract

The legrest can extend within a range of 7 inches.



● Legrest down

Press this button on the remote control to make the legrest go down. Release the button to lock in position.



● Arm massage

Put your arms between the arm airbags before they are filled with air.



● Legrest up and backrest down

Press this button to make the legrest go up and the backrest go down. Release to lock in position.



● Shoulder massage

Put your shoulders between the shoulder airbags before they are filled with air.




● Backrest up and legrest down

Press this button to make the backrest go up and the legrest go down. Release to lock in position.



● Leg massage

Multiple airbags will provide massage to the legs to promote blood circulation.



● ZeroG mode

Press this button to go into ZeroG mode. The backrest will go down and the legrest will go up. Release to lock in position.



● Foot massage

The foot rollers will provide massage to the feet, which helps the foot reflex zones.



- Warning**
1. Ensure no people, pets, or objects are within range of the chair.
 2. When adjusting the backrest, be careful of the gap between the backrest and armrest.
 3. When the legrest is moving, do not stand on it or leave the seat.




Massage operation

- 1 Make sure there are no people, pets, or objects within range of the backrest and legrest.
- 2 Return the backrest and legrest to the storage position by pressing the ON/OFF button.
- 3 After the chair resets, put the remote control back in its holder.
- 4 Put the pillow pads back on the backrest.



- 5 Switch off the power.
- 6 Unplug the power cord from the wall socket.
- 7 Clean the chair if necessary.



Attention!

- After using the chair, press the ON/OFF button to reset the chair.
- Unplug the power cord to avoid injuries.
- Keep the chair away from moisture and dirt.

Maintenance

Leather care

- **General cleaning**
Wipe chair with soft, clean dry cloth.
Do not use chemical products.
- **Leather care**
 - 1 Dip a soft cloth into 3-5% neutral cleaning solution. Squeeze the excess solution, then gently clean the surface.
 - 2 Wipe off the solution with a damp cloth.
 - 3 Wipe with dry cloth.
 - 4 Air dry. (Do not use air blower!)



Plastic parts

- 1 Dip cloth into cleaning solution, squeeze the excess and clean.
- 2 Wipe the chair again with damp cloth.
- 3 Air dry.



Remote Control

- 1 Wipe with a clean damp cloth. Do not use chemicals or rubbing alcohol.
- 2 Air dry.



Cloth material

- 1 Dip a soft cloth into weak neutral cleaning solution, then squeeze the excess and clean the surface.
- 2 Brush dirty spots with weak neutral cleaning solution.
- 3 Wipe the chair again with a damp cloth.
- 4 Air dry. (Do not use air blower!)



Before cleaning

- Unplug the chair from the wall socket.
- Make sure your hands are dry when plugging or unplugging the chair.



Problem	Possible reason	Troubleshooting
The chair doesn't operate after being turned on	The power switch may not be turned on.	Press down on the power switch.
	The massage program may not be selected.	Select the desired massage program.
	The cord may be damaged.	Contact your service agent.
	The cord may not be connected properly.	Ensure a proper connection.
	The fuse may have burned out.	Replace fuse of same specification.
	The inner circuit may have an issue.	Contact your service agent.
Continual beeping sound	The sound is coming from the air pump, motor, and mechanical parts.	This is normal, no action needed.
The left and right roller are different heights	The massage rollers work on shifts.	This is normal, no action needed.
The chair sounds much louder than normal	The massage parts may need a break due to long continual use.	Turn off the power for a half hour to let the chair rest.
	The inner parts may need to be serviced due to age. This is normal.	Contact your service agent.
The chair has powered down	The chair has powered down.	Connect the electricity again.
	The massage time is up.	Turn off the power for a half hour to let the chair rest.
Neither the backrest nor the legrest move	Check for any obstacles in the way. For safety reasons, sometimes the chair will shut itself off.	Remove any obstacles. Turn off the power for a half hour to let the chair rest.
The chair doesn't return to storage position	The chair has probably been in use for too long.	Contact your service agent.
The cord is hot	It may be overloaded.	Turn off the power for a half hour to let the chair rest.
	Something else may be wrong.	Contact your service agent.

Other

Other



Note: If the problem is not caused by the reasons stated above, turn off the chair and unplug it from the wall. Then contact your service agent.

Name	Massage Chair	Power	200W
Model	WS-19	NW/GW	Main body: 161/220 lbs Armrest: 31/44 lbs Legrest & footrest: 47/54 lbs
Input voltage	<input type="checkbox"/> AC120V, 60Hz <input type="checkbox"/> AC220V, 50/60Hz <input type="checkbox"/> AC240V, 50Hz	Packing size	Main body: 57.9 x 29.3 x 37 in Armrest: 44.5 x 16.5 x 24.8 in Legrest & footrest: 21.6 x 20.5 x 22.2 in (L×W×H)
Safety design	Class I		

